

DENTAL INSTRUMENTS CATALOG



www.bispnsintl.com
bispnsintl@gmtill.com
 +92-921 325 9595



table of contents



Diagnostics

1-6



Implantology

49-56



Oral Surgery

7-40



Micro Surgery

57




Oral & Maxillofacial

41-47



Endodontic

58-74



Luxatip

48



Periodontal

75-88



Prosthodontic

89-99



Accessories

128-129



Orthodontic

100-115



Instrument Management

116-127



International

www.bisonsintl.com
info@bisonsintl.com
bisonsintl@gmail.com

Diagnostics

Handle Options
Mirror Handles
Explorers
Micro Surgical Mirrors
Tweezers



Diagnostics

HOW TO USE THIS SECTION

Instrument Name & Pattern	23 Shepherd's Hook
Part Code of pictured instrument	EXS6XL
Available handle designs	Handle options : #4, #3, #2, #1

HANDLE OPTION

Most popular instruments available in all handle options.

DOUBLE END



SINGLE END




NOTE :
This is not a mouth mirror handle.
It's just only for the selection of instrument
handle which you want to use.


2.

Diagnostics

MIRROR HANDLE



MIRROR HANDLE SIGNATURE SERIES # GRIPS






DID YOU KNOW ?
Colour coded instruments rings are silicon
rubber rings for identifying and organising
instruments. Useful for maintaining proper
instruments order in cassettes and
identifying instruments by procedure.

3.

Diagnostics

PROBES

A double-ended instrument that combines two critical instruments required for every effective diagnostic set-up with an explorer on one end and a probe on the other, an Probe is excellent for dental examinations and charting.






4.

Diagnostics

MICRO SURGICAL MIRRORS


Double-sided-mirrors are manufactured from highly polished stainless steel and come in round-end or rectangular shapes.



DID YOU KNOW ?

Back surface
These mirrors have their reflective coatings on their rear surfaces.

Front surface
These mirrors have their reflective coating on their front surface, yielding undistorted imaging with no disturbing double images.



5.

Diagnostics

TWEEZER

Used to grasp or transfer materials in and out of the oral cavity. All dressing pliers are made of heavy-gauge, heat-treated stainless steel to prevent tip bending and misalignment.

London College
DPU17



Mariam
DP318



London College Premium
DP2



Mariam Premium
DP3



Locking Premium
DP18L
Locking Serrated
DP18LS



GP Holding Premium
EPL1
GP Holding
EPL1S



Brush Holding
BHT5



Perry
DP1



Curved Tweezer
DP90



Oral Surgery


Extraction Forceps

Root Elevators

Scissors & Needle Holders


Retractors

Periosteal Elevators




EXTRACTION FORCEPS

Most popular instruments available in all handle options.



S Standard
P Premium
E Ergonomic




1# Upper Anteriors
FX1S
FX1P
FX1E

7# Upper Premolars
FX7S
FX7P
FX7E

17# Upper Molars Right
FX17S
FX17P
FX17E

18# Upper Molars Left
FX18S
FX18P
FX18E

18A# Upper Molars
FX18AS



29# Upper Roots
FX29S
FX29P

29N# Upper Roots Narrow
FX29NS
FX29NP

30# Upper Roots
FX30S
FX30P

55# Separating Upper Molars
FX55S


56# Separating Lower Molars
FX56S

● Upper Anteriors ● Lower & Upper Molars ● Lower Anteriors ● Lower & Upper Premolars ● Lower & Upper Roots

EXTRACTION FORCEPS



S Standard
P Premium
E Ergonomic




41# Upper Roots
FX41S
FX41P

44# Upper Roots
FX44S
FX44P

51# Upper Roots
FX51S

51A# Upper Roots
FX51AS
FX51AP
FX51AE

67A# Upper Third Molars
FX67AS
FX67AP
FX67AE



76N# Upper Roots
FX76NS

89# Upper Molars Right
FX89S
FX89P
FX89E

90# Upper Molars Left
FX90S
FX90P
FX90E

13# Lower Premolars
FX13S
FX13P
FX13E

22# Lower Molars
FX22S
FX22P

● Upper Anteriors ● Lower & Upper Molars ● Lower Anteriors ● Lower & Upper Premolars ● Lower & Upper Roots

EXTRACTION FORCEPS



31# Lower Roots
FX31S



33# Lower Roots
FX33S
FX33P



33A# Lower Roots
FX33AS
FX33AP



33L# Lower Roots
FX33LS
FX33LP



73# Lower Molars
FX73S
FX73P
FX73E



74# Lower Anteriors
and Roots
FX74S



74N# Lower Roots
FX74NS
FX74NP
FX74NE




75# Lower Premolars
FX75S



79# Lower Third Molars
FX79S
FX79P
FX79E



45# Lower Roots
FX45S
FX45P



86# Lower Molars
FX86S
FX86P
FX86E



86A# Lower Molars
FX86AS
FX86AP



86B# Lower Molars
FX86BS
FX86BP



86C# Lower Molars
FX86CS
FX86CP




87# Lower Molars
FX87S




● Upper Anteriors ● Lower & Upper Molars ● Lower Anteriors ● Lower & Upper Premolars ● Lower & Upper Roots


EXTRACTION FORCEPS AMERICAN PATTERN



53R# Molars
First & Second Upper
Molar Teeth Right
F53R



53L# Molars
First & Second Upper
Molar Teeth Left
F53L



88R# Nevius
First & Second Upper
Molar Teeth Right
F88R



88L# Nevius
First & Second Upper
Molar Teeth Left
F88L




American Pattern




23# Cowhorn
First & Second
Lower Molars
F23



150# Cryer Universal
Upper Incisors, Canines,
Premolars & Roots
F150



151# Cryer Universal
Lower Incisors, Canines,
Premolars & Roots
F151



69# Tomes
Upper & Lower
Fragments & Roots
F69



150S# Peedo Universal
Upper Primary Teeth
& Roots
F150S



151S# Peedo Universal
Lower Primary Teeth
& Roots
F151S



Peedo Universal

DID YOU KNOW ?
All instruments should be cleaned and
sterilised in an open position.

● Upper Anteriors ● Lower & Upper Molars ● Lower Anteriors ● Lower & Upper Premolars ● Lower & Upper Roots


Universal Forceps



221# Universal
for Upper Roots
FX221



223# Universal
for Lower Roots
FX223



EXTRACTION FORCEPS # SECURE

Secure-Line Extraction Forceps

These forceps will allow you to complete any extraction and work more efficiently even if complications arise. If teeth break during extraction, these are the tools you need to rapidly, safely solve your problems. These forceps have special, allpurpose, jaws designed to cover a broad range of applications and have sharp pointed tips that allow separating the tips of roots from surrounding tissues by lightly twisting them and then securely trapping the extracted roots in the gaps between their jaws.




845.00# Lower Roots
SFX845.00



945.00# Lower Roots
SFX945.00



959.01# Lower Roots
SFX959.01




974.00# Lower Roots
SFX974.00



849.00# Upper Roots
SFX849.00




944.01# Upper Roots
SFX944.01




944.00# Upper Roots
SFX944.00



149.11# Upper Roots
SFX149.11




897.00# Upper Roots
SFX897.00



197.11# Upper Roots
SFX197.11



951.00# Upper Roots
SFX951.00



● Upper Anteriors ● Lower & Upper Molars ● Lower Anteriors ● Lower & Upper Premolars ● Lower & Upper Roots

EXTRACTION FORCEPS # ATRAUMATIC

Forceps "Blade Beaks"

- "Blade Beaks" smoothly penetrate in contact with the root.
- Large contact surface for a secure grip.
- Take care of the tissues and assure tooth fracture risk reduction.
- Better control over the root during the extraction operations.



1# Upper Centrals
& Canines
AFX1




7# Upper Premolars
AFX7



13# Lower Premolars
AFX13



51A# Upper Roots
AFX51A



74N# Lower Roots
AFX74N



79# Lower Third Molars
AFX79



Albrecht Extraction Forceps



1# Albrecht Premolar
FXA1

2# Albrecht Molar
FXA2

● Upper Anteriors ● Lower & Upper Molars ● Lower Anteriors ● Lower & Upper Premolars ● Lower & Upper Roots

Special Wisdom Router Extraction Forceps



22 1/2L#
Lower Molars & Wisdom Left
FX22 1/2LS
FX22 1/2LP



22 1/2R#
Lower Molars & Wisdom Right
FX22 1/2RS
FX22 1/2RP




Their offset design will allow you to reach around obstacles and protect the corners of patient's mouth whenever their mouth openings are restricted or their jaws have been firmly clamped in position. Their slim compact design allows readily accessing molar areas.

Oral Surgery

EXTRACTION FORCEPS BLACK LINE

The surface shows excellent sliding features to reliably counteract contamination. An incrustation of proteins etc. is not possible. The matt-finished black surface not only prevents reflection of the OT light, but also guarantees an optimal contrast in the OR environment during extraction.



● Upper Anteriors ● Lower & Upper Molars ● Lower Anteriors ● Lower & Upper Premolars ● Lower & Upper Roots



14.

Oral Surgery

EXTRACTION FORCEPS # PEEDO

STANDARD



5# Upper Incisors & Canines
FX5CS



1# Lower Incisors
FX1CS



2# Lower Primary Molars
FX2CS



3# Lower Roots
FX3CS



4# Upper Primary Molars
FX4CS



6# Upper Molars
FX6CS



7# Upper Roots
FX7CS



PREMIUM



7# Upper Roots
FX7CP



5# Upper Incisors & Canines
FX5CP



4# Upper Primary Molars
FX4CP



6# Upper Molars
FX6CP



1# Lower Incisors
FX1CP



2# Lower Primary Molars
FX2CP



3# Lower Roots
FX3CP




● Upper Anteriors ● Lower & Upper Molars ● Lower Anteriors ● Lower & Upper Premolars ● Lower & Upper Roots

15.

Oral Surgery

PROXITIP - LUXATING HYBRIDS

Designed with thin, precisely angled tips and a large diameter handle, the Luxating Hybrids provide improved apical and interproximal access for enhanced oral adaptation, aiding in the preservation of the buccal plate by severing the periodontal ligaments during atraumatic extractions.



Small Straight
EPTSMS
For front teeth
& upper jaw



Small Curved
EPTSMC
For lower jaw

Suggested Pair



Large Distal
EPTLGD
For large
distal luxating



Large Mesial
EPTLGM
For large
mesial luxating

Suggested Pair



Bernard Right
EPTBER6
For lingual /
buccal luxating



Bernard Left
EPTBER7
For lingual /
buccal luxating

Suggested Pair



Small Distal
EPTSDM
For distal luxating



Small Mesial
EPTSMm
For mesial luxating

Suggested Pair



Warwick James Left
EPTWJL



Warwick James Right
EPTWJR



Straight Spade
EPTSSP
For initial loosening
the desmodontal
fibres

Oral Surgery

ROOT ELEVATORS



Standard



Premium

Seldin



E34S
Thinner Tip & Narrow
than 34 Shank



E34



E2S
Lower Third Molar



E3S
Upper & Lower
Anterior Roots



E301W
For displacement of Deeply Seated Roots



E304W

Seldin



E1L
Upper & Lower Molar Roots



E1R



E4L
Smaller version of E1L & E1R
also useful for Impacted Third Molar



E4R



E1
Upper & Lower Anterior Teeth
and Roots



E2

Heidbrink Straight

Apexo



E301
For displacement of Deeply Seated Roots



E302



E71



E72
For displacement of Broken of Deep Seated Roots Tips



E73






E74

Miller Apexo

HANDLE OPTIONS



Oral Surgery ROOT ELEVATORS




18.



Oral Surgery

ROOT ELEVATORS

HANDLE OPTIONS



Also available in pouch set of 6 pcs.
LLSP6



19.

HANDLE OPTIONS



ROOT ELEVATORS

Oral Surgery


Oral Surgery

BONE RONGEURS

Manufactured from the finest grade materials for optimal performance and longevity, GDC features and extensive portfolio of rongeurs, including full round end cutting, end cutting and side cutting designs, engineered for trimming and recontouring alveolar bone and gross tissue removal.




Root Elevators Premium with Cassette



REPWC9

- Flohr #E1FP
- Flohr Right #E3FP
- Flohr Left #E2FP
- Cryer Right #E32P
- Cryer Left #E31P
- London Hospital #E20HP
- Coupland #EA1CP
- Winter Cryer Right #EWC12RP
- Winter Cryer Left #EWC12LP
- Cassette #IMSCDRE



Micro Bone Rongeurs


Small beaks for use in limited access areas and for trimming fibrous tissue or small bony spicules.



Oral Surgery

SCISSORS

GDC scissors are fabricated from special high-tempered surgical grade stainless steel providing longer lasting cutting edges and corrosion resistance. Designed to exacting specifications for balanced and smooth control. GDC offers an extensive selection including microsurgical, suture, crown and gold, bandage and general surgery scissors.



Iris Straight 11.5cm
S17



Iris Curved 11.5cm
S18



Micro Iris Straight 11.5cm
S26



Micro Iris Curved 11.5cm
S26C



Iris Side Curved 11.5cm
S27



Iris Angular 11.5cm
S7



Quinby Straight 12.5cm
S8S



Quinby Curved 12.5cm
S8



Goldman Fox Straight 13.0cm
S16



Goldman Fox Curved 13.0cm
S16C



Goldman Fox Double Curved 13.0cm
S10



Legrange Double Curved 12.0cm
S14

Oral Surgery

SCISSORS



Kelly Straight 16.0cm
S2



Kelly Curved 16.0cm
S1



Locklin Straight Handle 16.0cm
S12



Locklin Curved Handle 16.0cm
S11



Micro Heath for Suture Cutting 10.0cm
S25S



Heath for Suture Cutting 15.5cm
S25



Metzenbaum Straight 14.5cm
S28



Metzenbaum Curved 14.5cm
S29



Mayo Straight 14.5cm
S4



Mayo Curved 14.5cm
S3



Crown & Band Straight 12.0cm
SCGS



Crown & Band Curved 12.0cm
SCGC

Oral Surgery

SCISSORS

For Suture Cutting



Spencer 9.0cm
S13S



Spencer Side Curved 9.0cm
S13SA



Spencer 13.0 cm
S13



Spencer Side Curved 13.0cm
S13A

Marilyn Scissors



Straight 10.0cm
S3547



Curved 10.0cm
S3548



Straight 9.5cm
S3541



Curved 9.5cm
S3542



Dean Angular 16.5cm
S9



Noyes Angular 11.0cm
S42



Westcott Curved 11.0cm
S35



24.

Oral Surgery

SCISSORS



Castroviejo Straight 12.0cm
S32



Castroviejo Curved 12.0cm
S33



Castroviejo Straight 14.0cm
S39



Castroviejo Curved 14.0cm
S38



Noyes Straight 11.0cm
S30



Noyes Curved 11.0cm
S31



Micro Castroviejo Curved 18.0cm
SPV



Curved

Straight

25.

Oral Surgery

SCISSORS

Featuring tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.



Spencer TC Angular 10.0cm S3481



Spencer TC Angular 13.0cm S3482



Spencer TC Straight 10.0cm S3483



Spencer TC Straight 13.0cm S3484



Iris TC Straight 11.5cm S5082



Iris TC Curved 11.5cm S5083



Goldman Fox TC Straight 13.0cm S5080



Goldman Fox TC Curved 13.0cm S5081



Metzenbaum TC Straight 12.0cm S5054



Metzenbaum TC Cruved 12.0cm S5055



Crown & Band TC Straight 12.0cm S5038



Crown & Band TC Curved 12.0cm S5039

Oral Surgery

HEMOSTATS

Each hemostat is precision engineered from surgical grade stainless steel to exacting specifications featuring a secure ratchet and various working end geometries. Designed to optimise clinical efficiencies as a multipurpose instrument, the GDC hemostats are used for clamping off blood vessels, removing small root tips and grasping loose objects.



Kelly Straight 14.0cm H1



Kelly Curved 14.0cm H2



Crile Straight 16.0cm HRPS



Crile Curved 16.0cm HRPC



Mosquito Straight 12.0cm H4



Mosquito Curved 12.0cm H3



Mosquito Straight 1x2 12.5cm H5



Mosquito Curved 1x2 12.5cm H6



Baby Mosquito Straight 10.0cm HHS



Baby Mosquito Curved 10.0cm HHC



Allis Baby 12.5cm TFB



Allis 15.0cm TFS

Allis Forcep

Oral Surgery

NEEDLE HOLDER

Used to grasp and guide the needle when suturing



Mayo-Hegar Straight 16.0cm
NHMH



Mayo-Hegar Curved 16.0cm
NHMHC



Crile-Wood 15.0cm
NHCW



Adson Straight 17.0cm
NHA



Needle Holder 18.0cm
NHMH18

(Also available in 20.0cm #NHMH20)



Olsen-Hegar Needle Holder
with Scissors 17.0cm
NHOH



Derf Straight 12.5cm
NHD



Derf Curved 12.5cm
NHDC

Backhaus Towel Clamp



Backhaus 9.0cm
TC3



Backhaus 11.0cm
TC4



Backhaus 13.0cm
TC5

Castroviejo Needle Holder



14.0cm
Curved
NHCVC



14.0cm
Straight
NHCVS

28.

Oral Surgery

NEEDLE HOLDER TUNGSTEN CARBIDE

Needle holder feature tungsten carbide insert with a ground pyramidal profile to assure firm, positive needle grip. Rotation of needle and/or slipping of suture material is virtually eliminated



Mayo-Hegar 16.0cm
NH5042



Crile-Wood 15.0cm
NH5038



Halsey 13.0cm
NH5037



Derf 11.5cm
NH5032

Micro Castroviejo Needle Holder Tungsten Carbide Insert



Castroviejo TC Straight 14.0cm
NH5020



Castroviejo TC Curved 14.0cm
NH5021



Castroviejo TC Straight 18.0cm
NH5024R



Cheate Forcep 27.0cm
CF

Dressing Forceps



Gross Maier Straight 20.0cm
SDFS



Gross Maier Curved 20.0cm
SDFC



Foerster Straight 20.0cm
SDFF




Rogge 20.0cm / 30.0cm
RF20
RF30


29.

TISSUE FORCEPS

Used to grasp and stabilize soft tissue flaps during suturing.



MICRO TISSUE FORCEPS



DENTAL SYRINGE

GDC syringes are available in two styles, the cooke-waite and with large diameter finger rings for enhanced control, large window for enhanced carpule visibility and sharp harpoons for easy cartridge engagement.

Dental Syringe with colour coded ring



SYRACC1.8R

(Also available in different colours of ring)



Self Aspirating Syringe



SYRSA1.8

Intraligamental Syringe - Inject



SYRI1

"Suitable for upper and lowers, just turning the front barrel"



Not Aspirating

SYRI1.8
SYRI2.2

Aspirating

SYRA1.8
SYRA2.2

Aspirating

SYRF1.8

(Also available in European & American thread needle type)

ADVANTAGES OF INTRALIGAMENTAL ANESTHESIA :

- Automatic limitation of the power pressure.
- Anesthetizing of a single tooth.
- No numbness of the tongue, lips or cheek.
- Simultaneous treatment of teeth in different quadrants.
- Supports regional block anesthesia.
- Eliminates the risk of haematoma.
- Lower emotional impact on patient.


WHY INTRALIGAMENTAL IS BETTER

Because it's chrome free.
Because you can get a better sterilization.
Because you can use for long time
Because has anti-reflection surface.

OPERATING FEATURES:

- Fully autoclavable.
- Removable headpiece with angled threaded nozzle for the injection needle.
- Viewing chamber with safety sleeve for carpule cartridge.
- Handle servo dosing lever with optimal support for index finger for a gradual advance of the dosing plunger resetting key.
- No waiting time between injections and beginning of treatment.
- Excellent depth of anesthesia combined with shorter duration.
- Much less pain than standard injections (infiltration anesthesia in hard palate, intrapulpal or intraseptal injections or for papillectomy).
- A smaller quantity of given anaesthetic (~0,06 ml for single injection/click.)
- Needle suggested: 30G.

SUCTION TUBES & ADAPTOR



- Coupland Tip Size 1.5mm ASPCT1
- Coupland Tip Size 2.0mm ASPCT2
- Coupland Tip Size 2.5mm ASPCT3
- Coupland Tip Size 3.0mm ASPCT4
- Coupland Tip Size 5.0mm ASPCT5

- Frazier Tip Size 1.0mm ASPFR4
- Frazier Tip Size 1.5mm ASPFR6
- Frazier Tip Size 2.0mm ASPFR8
- Frazier Tip Size 2.5mm ASPFR10
- Frazier Tip Size 3.0mm ASPFR12
- Frazier Tip Size 3.5mm ASPFR14

To connect aspirator with light vacuum

Aspirator Adaptor with Silicon Tubing PTAA

BONE CHISEL

Used to scoop bone and make trough



Gardner 3mm
CG1

Gardner 4mm
CG2

Chandler 5mm
CC7

Chandler 7mm
CC9

MALLETS

Used with a chisel to split teeth and remove bone




Mead 2 17.0cm
MAL2

Mead 1 19.0cm
MAL1

BONE FILE

Bone File used for final smoothing of bone



Miller
BF52S

Miller 52
BF52

Howard 12
BF12

Miller 64
BF64

Oral Surgery

SCALPEL HANDLES

GDC scalpel handles feature balanced designs and are manufactured from the finest grade surgical stainless steel for optimal control and enhanced longevity.



Scalpel Handle Straight 14.5cm
10-130-5EM



Scalpel Handle Angular 14.5cm
10-130-5AE



Scalpel Handle Straight 14.5cm
10-130-5EM #4



Scalpel Handle Angular 14.5cm
10-130-5AE #4



Bard Parker No. 3 12.5cm
10-130-03




Bard Parker Style with
Metric Ruler No. 3 12.5cm
10-130-03E



No. 4 13.5cm
10-100-04E



No. 7 16.5cm
10-130-07E



Moon's Probe
P13M



Orientable Scalpel Handle
10-130-180

Double Side Blade
Scalpel Handle
10-130-5DE

34.

Oral Surgery

CHEEK RETRACTORS / TISSUE RETRACTORS

Used to hold mucoperiosteal flaps, cheeks, lips, and tongue away from the surgical area.

Austin Cheek Retractors



Regular
TRA4

Fork
TRA2



Tissue Retractor
with Handle
TRA1

Halfmoon
TRA3



Farabeuf
12.0cm
CRF

Mouth Gags

Molt



Heister 13.0cm
HMG



Large 14.0cm
MGA



Small 11.0cm
MGC

Hartmann Retractor



15mm 72x24mm
CRH

Weider Cheek & Tongue Depressor



Large 15.0cm
CRW2

Small 14.0cm
CRW1

Orringer Retractors



60x115mm Large
CRO3

Also Available in CRO2 Medium & CRO1 Small

35.

CHEEK RETRACTORS / TISSUE RETRACTORS




Bishop 13.5cm
CRB

Vestibulum 14.0cm
CRV

Martin Simplex 10.5cm
CRMS

MSY 11.5cm
MSY



Columbia 14.0cm
CRC

Sternberg 16.0cm
CRCS

Bruenings 19.0cm
TRB

University of Minnesota
14.0cm
CRIM

Cheek Retractor with Suction




Fork
TRA2S

Half Moon
TRA3S

With Suction Large 15.0cm
CRW2S

University of Minnesota
with Suction 14.0cm
CRMS

CHEEK RETRACTORS / TISSUE RETRACTORS



30 x 10
CRL1

40 x 10
CRL2

40 x 13
CRL3

50 x 10
CRL4

14 x 17
CRMD1

20 x 22
CRMD2


26 x 30
CRMD3

35 x 11
CRR

Langenbeck 21.5cm

Middledorf 19.0cm

Ramus 21.5cm



56 x 11
CRK2

40 x 11
CRK1

55 x 11
CRK4

40 x 11
CRK3

Kocher 21.5cm

Kocher 21.0cm

CHEEK RETRACTORS / TISSUE RETRACTORS



Ragnell - Davis

Senn - Miller

Appuntito



Flexible Con Molla 16.5cm

Gillies Retractor

PERIOSTEAL ELEVATOR

For reflecting and retracting the mucoperiosteum after incisions of the gingival tissue.



Hirschfeld
P20
#6

Glickman
P24G
#6

Kramer-Nevins
PKN1
#6



Prichard
PPR3
#6

Benque
PBQ
#6

Buser
PPBUSER
#6



Seldin
P23



Goldman Fox
P14
#6

P149
#6

P8D
#6

PH2
#6



Howarth
PH

PERIOSTEAL ELEVATOR



PPLABANCA2
#6

P856
#6

Molt
P883
#6

PHP3
#6

Freer
P884
#6



Molt-9 Large
P9
#3

Molt-9 Small
P9S
#3

Freer
P16
#3

PPLABANCA2
#3

BEAL1
#3

Kramer Naviens
PKN
#3


Oral Maxillofacial

Retractors

Trauma Kit Instruments

Cleft Lip & Palate Instruments

Impaction Kit Instruments



Oral Maxillofacial

BONE HOLDING, RETRACTORS & OTHER INSTRUMENTS

Asch Nasal
Septum
Forcep



Large
ANS

Walsham
Nasal Septum
Forcep



Left
WNSL

Right
WNSR

Rowe Maxillary
Disimpaction
Forcep



Left
RMDL

Right
RMDR

Nasal Chisel
with Guard



Left
NCGL

Right
NCGR

Ferguson
Mouth Gag



Large
FMGA

Small
FMGC

Hayton William
Forcep



Left
HWL

Right
HWR

42.

Oral Maxillofacial

BONE HOLDING, RETRACTORS & OTHER INSTRUMENTS



Bone Holding Forcep (Mandible)
BHM



Channel Retractor Regular
CHR



Kilner Retractor C Shape
KRC



Chin Retractor
CR



Pterigoidous Chisel
PC



Nasal Septum
Osteotome with
Guard U Shape
NSOU



Bristow Elevator
BE



Bone Elevator
Hohman Type
BEH



Micro Miller Bone
File Straight
BFMS



Barsky Pharyngeal
Flap Elevator
BPF

43.

Oral Maxillofacial

BASIC BONE PLATING INSTRUMENTS



Modelling Lever
ML



Plier
BPP



Bone Cutting Forcep Single Action
RBC



Zygomatic Bone Awl
ZBA

Suggested Trio



Suggested Trio



Wire Twister
WTL



Bone Holding Forcep Pointed
BHP

Ordinary Screw Driver
OSD1.5
OSD2.0
OSD2.5

Self Holding Screw Driver
SSH1.5
SSH2.0
SSH2.5




Trauma Kit
TRAUMA

Suggested Trio



SS Drill Bite for Micromotor
DSS1.5
DSS2.0
DSS2.5



Orbital 4 Holes P4HO1.5SS P4HO2.0SS
Left 4 Holes P4HLL1.5SS P4HLL2.0SS
Right 4 Holes P4HLR1.5SS P4HLR2.0SS
16 Holes P16H1.5SS P16H2.5SS

44.

Oral Maxillofacial

CLEFT LIP & PALATE INSTRUMENTS



Palate Dissector Right
CPDR



Palate Dissector Left
CPDL



Cleft Palate
Raspatory Up
CPRU



Cleft Palate
Raspatory Down
CPRD



Senn - Miller Retractor (Kilner)
TRSM1



McIndoe Cleft Palate Raspatory
MCPR



Cleft Palate Elevator
CPE2



Spoon Shape Elevator
SSE



Cleft Palate Elevator
CPE1




Miller Bone File Curved
BF52

45.

Oral Maxillofacial


CLEFT LIP & PALATE INSTRUMENTS




46.

Oral Maxillofacial

IMPACTION KIT



47.



LUXATIP

GDC Luxatip Luxation Instrument

Over recent years, there has been an increased demand for reduced trauma during extraction procedures. The use of elevators has become widespread for the loosening of teeth before the application of forceps. GDC Luxation instruments tips are more finely tapered than elevators tips. Their characteristics are a combination of a Perioste and an elevator. The sharp tip enables more efficient cutting of the periodontal ligament, whilst the tapering shaft assists with controlled expansion of the bony socket. The extraction operation may be performed with the minimum amount of tissue damage. Healing is therefore quicker and the operation is less distressing for the patient and surgical team alike.

The GDC Luxatip Kit comprises of seven Luxation instruments, a straight, curved, inverted curved & contra blade version in Three different widths.

Also available in dual edge blades

- L2S Straight Blade- Green- 2 mm- Apical and/or Interproximal
- L3C Curved Blade- Black- 3 mm General use/Starter
- L3S Straight Blade- Grey- 3 mm Interproximal
- L3CA Contra Angle- Blue- 3 mm Lingual and/or distal-Molars
- L3IC Inverted Curved- Light Blue- 3 mm- General, Lingual and/or Distal
- L5C Curved Blade- Brown- 5 mm- Large molar roots-General Use
- L5S Straight Blade- Yellow- 5 mm- Interproximal
- DE3 dual edge/straight grey 3 mm/1.5 mm
- DES dual edge/straight brown 5 mm/3 mm

Deep fractured/Decayed Roots



Implantology



Osteotomes

Sinus Lift Instruments

Caliper & Chisel

Bone Injector & Collector



Bone Mill & Crusher



OSTEOTOMES


Osteotomes Bone Pushers Straight

Allows for the transport of bone graft chips from the implant cavity or for condensing them apically.



Osteotomes Bone Pushers Angulated



Designed to condense bone material and/or autogenous bone under the sinus. Minimise the risk of perforation of the mucous membrane in sinus elevation procedures.





OSTEOTOMES CONVEX & CONCAVE

Osteotomes Tapered Convex Tip Angulated

Specially designed for closed sinus elevation operations in the lateral upper jaw region or in difficult to access maxillary areas.



Osteotomes Tapered Concave Tip Angulated



Implantology

BONE MILL & SCRAPER

Bone Mill

For gaining autogenic bone material.

In order to correct bone defects, an intraoral harvesting of bone is often necessary. Pieces of bone may be gained with trephines or chisels. It is the aim of the Bone Mill Forceps to cut bone material into small pieces so that a maximum volume of bone chips can be achieved. The "grainy" consistency of bone chips created will ensure that it is adapted most favourably to osseous embedding.



Bone Crusher
IMPBMS



Bone Cracker for Large Bone Pieces
IMPBMR



Bone Mill
IMPBM

Bone Scraper

Collecting autologous bone material.

In oral surgery bone replacements and bone structures are often combined with autologous bone pieces during augmentations in order to use their osteoinductive effect. The Bone Scraper is the perfect instrument for an easy and quick collection of autologous bone structures. It enables you to scrape, collect and transplant the patient's own bone. Everything can be done without the usage of bone filter, trephine bur, saw or bone mill.

The blade makes a collection of cortical and spongiose bone possible. The bone is collected in a chamber during the scraping and at the same time the bone is being mixed with blood. The bone material can be implanted directly out of the Bone Scraper which has been sterilized before.

DID YOU KNOW ?

Bone Scraper Blades can be purchased separately.



IMPBSB2



IMPBSB1



IMPBS



IMPBC

Bone Scraper Blade

Bone Scraper

Implantology

IMPLANTOLOGY CALIPER & BONE COLLECTOR

Bone Caliper



BC15
15.5cm



BC35
15.5cm

Castroviejo Caliper



CLC40L
18.0cm



CLC40C
18.0cm



CLC20C
0.5cm



CLC20L
0.5cm

Bone Collector



BC1331/1F

BC1331/1
18.5cm



BC1331/2F

BC1331/2
18.5cm

Trephine Size (in mm) :
3 | 4 | 5 | 6 | 7 | 8 | 9 | 10




Trephine Kit
IMPTS8




Drill Size (in mm) :
2 | 2.5 | 2.8 | 3.2 | 3.65 | 4.2 | 4.8 | 5.2

Drills Kit
IMPDS8

TISSUE PUNCH



Instrument to measure the Depth of bony tunnel



PRF BOX
Palatelet Rich Fibrin
PRF

OSTEOTOMES BONE CHISEL



OCHSENBEIN CHISEL




BONE PUNCH



Ratchet Wrench Universal
RW

BONE INJECTORS



MALLETS

Used with a chisel to split teeth and/remove bone




Mead 1
MAL1
19.0cm

Mead 2
MAL2
17.0cm

SINUS LIFT INSTRUMENTS

Featuring a wide range of sinus lift instrument designs, GDC offers a variety of size, working end and angulation options to provide optimal access to separate and reflect the schneiderian membrane from the maxillary bone and to elevate the membrane during lateral at sinus lift procedures.



IMPSL1

IMPSL2

IMPSL3

IMPSL4

IMPSL5



IMPSL6

IMPSL7

IMPSL8

IMPSL9


MICRO SURGERY INSTRUMENTS



Castroviejo Needle
Holder 18.0cm
NH5024R



Castroviejo Scissor
18.0cm
SPV



Micro Surgery
Instruments With
Cassette
MSIWC6



Micro Tissue
Forcep Straight
TPSSTMBH



Micro Scalpel
Handle 12.0cm
1-015



Papilla Elevator
PH26M



Glickman
P24G

Endodontic

Composite Instruments

Spoon Excavator

Plugger / Condenser / Burnisher

Root Canal Plugger & Spreader

Rubber Dam Clamp, Forcep & Punch

Endodontic

COMPOSITE INSTRUMENTS

Plasma Black



CIPB1
#8

CIPB2
#8

CIPB3
#8

CIPB4
#8

CIPB5
#8

CIPB6
#8

Titanium Blue



CIB1
#8

CIB2
#8

CIB3
#8

CIB4
#8


CIB5
#8

CIB6
#8

Endodontic


COMPOSITE INSTRUMENTS

Titanium Gold



CIG1 #6 CIG2 #6 CIG3 #6 CIG4 #6 CIG5 #6 CIG6 #6

Titanium Blue




CIB1 #6 CIB2 #6 CIB3 #6 CIB4 #6 CIB5 #6 CIB6 #6

Endodontic

COMPOSITE INSTRUMENTS


Titanium Blue # Anterior



TNCIGFT3 #6 TNCIPCL #6 TNCIPCM #6 TNCIPCS #6 TNCIGFTMI1 #6

Also available in Plasma Black

Titanium Blue # Posterior




TNCFIS/M #6 TNCFIM/L #6 TNPLG5A #6 TNPLGOT #6 TNPLGH3 #6

Also available in Plasma Black

COMPOSITE INSTRUMENTS

Titanium Goldstein Flexi Thin



CIGFT1
#6

CIGFT2
#6

CIGFT3
#6

CIGFT4
#6

CIGFTMINI3
#6

CIGFTMINI1
#6

Anterior Kit

TNCIGFT3
Flexible, reversed, flared paddle design for shaping and placement of Class III and IV restorations.

TNCIPCL
For final placement in Class I and II restorations. The larger, round ball end is used for condensing and shaping in Class I and II restorations and on lingual surfaces of anterior teeth.

TNCIPCM
For small pits and fissures, as well as placement and condensing with limited access.

TNCIPCS
Mico-Mini for extremely small pits and fissures.

TNCIGFTMI1
Mini versions of the TNCIGFTM1 for small pits and fissures, tunnel preparations or minor tooth defects on lower anteriors.

Posterior Kit

TNCFIS/M
Rounded cone-shaped paired instrument designed to provide improved contact forming for small/medium Class II restorations.

TNCFIM/L
Rounded cone-shaped instrument to provide improved contact forming for medium/large Class II restorations.

TNPLGSA
Small, round, inverted-cone plugger for use with condensable composite material in posterior restorations.

TNPLGOT
Rhomboid-shaped plugger for use with condensable composite material in posterior restorations.

TNPLGH3
Rectangular-shaped plugger for use with condensable composite material in posterior restorations.

Gold Flexi Thin

CIGFT1
Small universal style with rounded plugger tip and a narrow paddle for initial placement and contouring of Class I, II and III restorations.

CIGFT2
Large universal style for final placement and contouring of Class I, II and III restorations.

CIGFT3
Flexible, reversed, flared paddle design for shaping and placement of Class III and IV restorations.

CIGFT4
Flexible, paired, offset, paddle-shaped bladed for placing and shaping material on posterior, mesial and distal surfaces. Reverse angle is also useful for placing and shaping anterior bonded restorations.

CIGFTMINI3
Mini version of the CIGFT3. Can also be used for packing gingival retraction cord.

CIGFTMINI1
Mini version of the CIGFT1 for small pits and fissures, tunnel preparations or minor tooth defects on lower anteriors.

SPOON EXCAVATOR

GDC endodontic excavators feature extended length shanks and are meticulously hand crafted from the highest grade stainless steel. They help the clinician access the pulp chamber for easy and efficient removal of dentin. The result is a superior clinical outcome and improved efficiency.



EXC125/126
#3 #6

EXC127/128
#3 #6

EXC129/130
#3 #6

EXC131/132
#3 #6

EXC133/134
#3 #6

EXC14
#3 #6

Oval Spoon



EXC17H
#3 #6

EXC18H
#3 #6

EXC19H
#3 #6

EXC17WH
#3 #6

EXC18WH
#3 #6

EXC19WH
#3 #6

Endodontic

SPOON EXCAVATOR



EXC31L #3 #6
EXC32L #3 #6
EXC33L #3 #6
EXC31W #6
EXC816 #3
EXC817 #3
EXC819 #3



EXC820 #3
EXC210/211 #3 #6
EXC212/213 #3
EXC153/154 #3 #6
EXC38/39 #6
EXCGL2 #6

Endodontic

PLUGGER / CONDENSERS & PLASTIC FILLING INSTRUMENTS

Used to pack filling materials into cavity preparations.



PLG1 #3 #4
PLG2 #3 #4
PLG3 #3 #4
PLG4 #3 #4
PLG1L #3 #4
PLG2L #3 #4
PLG3L #3 #4
PLG4L #3 #4

Plastic Filling Instruments



PFGFT1 #3
PFGFT2 #3
PF2 #3 #4
PF21 #3 #4
PF1 #3
PF12 #3
PF13 #3
PLGH3 #3 #4

BURNISHER & CARVERS

Designed to condense, smooth, carve and polish amalgam



BB21B #3 #4 BB21 #3 #4 BBL3 #3 PKT3 #3 #4 BB22/23 #3 #4 BB18 #3 #4 BB27/29 #3 #4 TB2 #3 #4

Carvers

Used to carve anatomical features and trim excess materials.



Discoid-cleoid CD89/92 #3 Discoid-cleoid CD4/5 #3 #4 Interproximal Off-angle CVIPCOA #3 Hollenback CVHL3S #3 #4 Wall W-Shaped CVHL3W #3 Ward's CVWR1S #3 #4 Diamond Shape Frahm CVFR2/3H #3 #4 Diamond Shape Frahm CVFR2/3 #3


MARGINAL TRIMMER

Used to produce proper bevel on enamel margins.



MT26H #6 MT27H #6 MT28H #3 #4 #6 MT29H #3 #4 #6


ENDO GUNS FOR RETROGRADS APPLICATIONS SYRINGES



MTA4865 13.0cm MTA4866 13.0cm MTA4869 13.0cm

PLACEMENT INSTRUMENTS

Used to place base or liner within cavity preparations.



Calcium hydroxide or glass ionomer base/liner placement instrument. Also useful as a small burnisher.

Calcium hydroxide or glass ionomer base/liner placement instrument combined with a mini-spatula for efficient mixing.

Calcium hydroxide or glass ionomer base/liner placement instrument.


PICH6 #3 SP60616 #3 PICH7 #3 PICH8

Endodontic

CHISELS AND HOE

Used to refine and cavity preparations. Forming lines angles on anterior preparations.

Enamel Chisel



CP7
#3 #4

CP7/8
#3

Wedelstaedt



CP1/2
#3

CP3/4
#3

CP5/6
#3

Hatchet



CP53/54
#3 #4 #6

Hoe



CP24
#3 #4

MATRIX RETAINER



Ivory
MRI



Tofflemire Universal
MRTU



Squiavand Narrow
MR4972S



Squiavand Wide
MR4972L

AMALGAM CARRIER

To carry and dispense amalgam filling materials.



AC504



Large/Regular
AC5202



Jumbo/Regular
AC5203



Mini/Regular
AC5201

Amalgam / Bone well
WA

Stainless steel, with removable non-slip base ring.
Designed for easy amalgam/bone graft handling.


Endodontic

ISO*SIZED PLUGGER COLOUR CODED


When performing vertical condensation, choose from ISO-sized anterior or posterior. These pluggers are handcrafted to exact specifications using the finest materials.

* DID YOU KNOW ?

ISO Pluggers have an ISO colour-coded ring on the handle for ease of identification during usage.
#30 Blue #40 Green #50 Yellow #60 Red



LUKS




Endodontic

SPREADERS, HEAT CARRIERS

ISO sized spreaders with color coded handles for easy identification, heat carriers, and a variety of other popular spreaders. Our spreaders help the clinician gain access into the most difficult canals for optimal compaction of filling material during lateral condensation.



RCSGP1 #3 #7
RCSGP2 #3 #7
RCSGP3 #3 #7
RCSMA57 #3 #7



RCSD11 #3 #7
RCSD11T #3 #7
RCSD11S #3 #7
RCSD11TS #3 #7

ROOT FRAGMENT STIEGLITZ

Post and Silver Point Removal forceps



RFSTR RF45 RF90

NERVE-CANAL PLIERS

Nerve-Canal Pliers for grasping fractured root-canal instruments or silver pins. The very fine, concave cutout on the working tip allows to remove even deep-lying fractured files.



RFES



RFEC


Endodontic

RUBBER DAM

Rubber Dam Clamps # Gold Line



Rubber Dam Clamps # Black Line



- POINTS OF PERFORMANCE
- Which resists corrosion while aiding in flexibility and memory
 - A dark coating allowing for improved intraoral contrast and visual acuity
 - Matt finish for decreased light reflection and glare



Endodontic

RUBBER DAM

Rubber Dam Clamps # Matt Finish



2A
RDC2A



3
RDC3



4
RDC4



12A
RDC12A



13A
RDC13A



14A
RDC14A



201
RDC201



202
RDC202



205
RDC205



210
RDC210



211
RDC211


Rubber Dam Clamps Set



Gold
RDCOBG11



Black
RDCOBB11



Matt Finish
RDCOB11

Rubber Dam Clamps Peedo # Matt Finish



A1
RDCA1



A2
RDCA2



A3
RDCA3



A4
RDCA4



A5
RDCA5



A6
RDCA6

DESCRIPTIONS

2A
Flat jaws for larger premolars.

3
Flat narrow jaws for small lower molars.

4
Flat narrow jaws for small upper molars.

12A
Serrated offset jaws for lower right molars. Paired with 13A

13A
Serrated offset jaws for lower left molars. Paired with 12A

14A
Deeply festooned jaws for partially erupted or irregularly shaped molars.

201
Slightly festooned jaws for upper molars.

202
Universal molars.

205
Slightly festooned jaws for large upper molars.

210
Universal double bowed anterior clamps.

211
Universal double bowed anterior clamps.



Winged Clamps



Wingless Clamps

Endodontic

RUBBER DAM CLAMP, FORCEP & PUNCH



Ivory 17.0cm
RDF



Stoke 17.5cm
RDF4E



Ainsworth 17.0cm
RDPA



Rubber Dam Punch 16.0cm
RDP




Matrix Ring Holding Forcep 17.0cm
MRDF




Dental Dam Kit Adult DDK
Dental Dam Kit Peedo DDKP

Endodontic

RUBBER DAM



ART Instruments



The State of the art - ART Instruments

ART restoration step-by-step using a high-viscosity glass-ionomer (Ketac Molar EasyMix). Dentine carious cavity in tooth LLE.

(a) Notice the discoloration around the cavity opening. It indicates that the carious lesion has extended under the enamel. This unsupported enamel is demineralized and will break off easily under slight pressure.

(b) Opening of cavity further for access with the Enamel Access Cutter (ART instrument).

(c) Carious tissue removal using a small and sharp excavator.

(d) The conditioner is applied to the cleaned cavity, pits and fissures with a cotton wool patted.

(e) Carefully dried cavity before placing the high-viscosity glass-ionomer.

(f) Firm finger pressure has been applied over the occlusal surface and excess filling material has been removed with a carver instrument. The ART restoration presented after the bite has been adjusted showing the sealant restoration: the cavity is restored, fissures are sealed.

Periodontal

Curettes


Scalers

Probes

Knives & Chisels

Periotomes





STANDARD

Gracey Curette Standard
Originally designed to be finishing currettes for fine scaling and root planing. The combination of this unique blade and 9 different shank designs for specific tooth surfaces provide improved adaptation and deposit removal in scaling procedures.



MINI FIVE

Gracey Curette Mini Five


- Terminal shank is 3mm longer than Standard Gracey
- Compared to the standard Gracey curette, blade is 50% shorter and 10% thinner
- Good for scaling in deep, narrow pockets

HANDLE OPTION




#3 #6 #7

GRACEY CURETTES # STANDARD



SG1/2 #3 #6 SG3/4 #3 #6 SG5/6 #3 #6 SG7/8 #3 #6 SG9/10 #3 #6 SG11/12 #3 #6 SG13/14 #3 #6 SG15/16 #3 #6 SG17/18 #3 #6


GRACEY CURETTES COLOUR CODED # STANDARD



SG1/2 #7 SG3/4 #7 SG5/6 #7 SG7/8 #7 SG9/10 #7 SG11/12 #7 SG13/14 #7

Also available Gracey Currettes Colour Coded # RIGID

GRACEY CURETTES COLOUR CODED # MINI FIVE



SAS1/2 #6 #7 SAS3/4 #6 #7 SAS5/6 #6 #7 SAS7/8 #6 #7 SAS11/12 #6 #7 SAS13/14 #6 #7 SAS15/16 #6 #7

Periodontal


CURETTES

Mc CALL CURETTES




SM11/12 #6 SM11/12S #6 SM11/12A #6 SM13/14 #6 SM13/14S #6 SM17/18 #6 SM19/20 #6

LANGER CURETTES



SL1/2 #6 SL3/4 #6 SL5/6 #6
Langer 1/2 for mandibular posterior.
Langer 3/4 for maxillary posterior.
Langer 5/6 for maxillary and mandibular anterior.

COLUMBIA CURETTES



Columbia 2R/2L anterior application
Columbia 4R/4L posterior application

78.

Periodontal

SCALERS

SCALERS



Townner-Jacquette SU15/30 #3 #4 #6 Springfield SSP2/3 #3 #6 Hoe Orban Mesial/Distal SO6/7 #3 #6 Hoe Orban Buccal/Lingual SO8/9 #3 #6 Posterior SN135 #3 #4 #6 Anterior SN137 #3 #4 #6 Crane-Kaplan SCK6 #6



Hygienist SH6/7 #6 Anterior SCNEVI1 #6 Posterior SCNEVI2 #6 Cumine SCM152 #3 #6 Posterior S204S #3 #6 Remington SRM3/4 #6 Surface MDY652 #3 #6

79.

Derby Perry
Chisel

MANIPAL SCALERS



Posterior
S129
#3 #6

Hoe Urban
SO7/9
#3

Bifid
SB4/8
#3 #4

Derby
SDP5/6K
#6

SM137


SM2/3

SM5/6

TUNNELLING KNIVES

FURCATION CURETTES

Designed for scaling root concavities and furcation.



Tunnelling Knife 1

Tunnelling Knife 2

TKN1
#6

TKN2
#6

The shank of the TKN1 has a smaller angle so it is suitable for tunnelling work in less complex areas.

The small radius shank of the TKN2 allows incision in convex areas with relatively large undercuts.

Quentin Buccal/Lingual

Quentin Mesial/Distal

SQBL2
#6

SQMD2
#6

SUB GINGIVAL SCALERS



SF1

SF2

SWV1


SWV2

SEC1

SEC2

SCU61

SUPRA GINGIVAL SCALER



SU15

SSP1

SB3

SB2

SCI1

SB6

SDP5

SDP6

PERIO PROBES

Used to assess periodontal pocket depths, attachment levels, anatomy configurations and gingival bleeding.



Thin William
POW6
#3 #4 #6


Goldman Fox
PGF6
#6

University of North Carolina
PCPUNC15
#3 #6

Marquis Colour Coded
PCP12
#6

Screening Probe
PCP11.5B
#3 #6

PERIO PROBES




EXP23/QW6
#6

EXP23/UNC15
#6

EXP23/11.5
#6

University of North
Carolina/William
PCPUNC15/W6
#3 #6

Naber's
Colour Coded
PQ2N
#3 #6 #7



Thin William/
Screening Probe
PCPQOW/11.5B
#3 #6

Goldman Fox/
William Probe
PCPGF/W6
#6

UNC15/
Screening Probe
PCP15/11.5B
#6

CP12/
Thin William
PCP12/QOW6
#6

CP12/
Screening Probe
PCP12/11.5B
#6

HIRSCHFELD FILE

Used to crush and remove heavy calculus deposits.

PERIODONTAL FILES

Used interproximally. Files surfaces on both sides allow for push or pull application.



Buccal/Lingual
FH3/7
#6

Mesial/Distal
FH5/11
#6

Anterior
FH9/10
#6

Sugarman
FS1/2S
#6

Sugarman
FS3/4S
#6

Schluger
FS9/10S
#6

For removal of granulomatous tissue in supra-bony and intra-bony pockets.

OCHSENBEIN CHISEL

Useful for reshaping bone. The semi-circular notches are beveled on each side. Can also be used to reflect flaps or remove secondary palatal flaps.



CO1

CO2

CO3
#6

CO4
#6

C36/37
#6

RHODES BACK ACTION

Designed for recountouring bone. Particularly valuable for removal of 'widows' peaks on difficult to reach distal surfaces of posterior teeth.

PERIODONTAL KNIVES

Goldman-fox



KGF8
#6

KGF9
#6

KGF11
#6

Buck



KB3/4
#6

KB5/6
#6

Crane Kaplan



KCK3
#6

Corn Suture Plier



SP20SLEM
18.0cm

SPM20
15.0cm

SP20

Orban



KO1/2
#6

KO12KP03A
#6
30% thinner

Kirkland



KK15/16
#6

Periodontal

Tissue Nipper



NIPS

Paquette
Blade Handle



KPQ

Goldman Fox
Pocket Marker



Left
PMGF1


Right
PMGF2

Surgical Curette Lucas



CL84 #3 #4 #6
CL85 #3 #4 #6
CL86 #3 #4 #6
CL87 #3 #4 #6
CL88 #3 #4 #6

Bone Curette Hemingway



CH876 #3 #4
CH877 #3 #4
CH878 #3 #4
CH879 #3 #4


Periodontal

PERIOTOMES # FLEXI

This product is used for minimizing damage to the alveolar bone as it helps cut the periodontal ligament to detach teeth from the periodontal ligament for implant immediately after tooth extraction or general tooth extraction.


CAUTIONS

PERIOTOME FLEXI is used for cutting off only the periodontal ligament. Do not lift it up or twist with excessive force such as ELEVATOR.



PTF1

PTF2




PTF3

PTF4

PTF6

PERIOTOMES NON FLEXIBLE




P1

P2


P3

P4



VISTA TUNNELING PROCEDURE

The first incision is made at the base of the sulcus, then sever the attachment from the osseous crest raise the mucoperiosteal pouch with a sub periosteal dissection and follow with a curette to separate the Vista from the osseous crest.



VISTA1 VISTA2 VISTA3 VISTA4 VISTA5 VISTA6

- Corn Suture Plier #SP20
- Scaler #S204S
- Surgical Curette Lucas #CL86
- Castroviejo Scissor Curved #S33
- Scalpel Handle Angular #10-130-5AE
- Perio Probe Thin William #POW6
- Suction Tubes Coupland #ASPCT3
- Periosteal Elevator Benque #PBQ
- Micro Tissue Forceps Gerald Straight Serrated #TPG1
- Micro Tissue Forceps Gerald Curved Serrated #TPG2
- Needle Holder Mayo-Hegar Straight #NHMH
- Tissue Forcep Adson 1x2 #TP42
- Instrument Management System #IMSCA20E



Perio Surgical Instruments with Cassette
PSIWC12

Prosthodontic

Wax & Modelling Carvers

Cement Spatulas


Wax Knives & Spatulas

Impression Trays


Crown Removers




WAX AND MODELLING CARVERS




PKT CARVERS COLOUR CODED



PKT CARVERS



WAX KNIVES & SPATULAS



HEIDEMANN SPATULA



PFIHS1
#1 #6

PFIHS2
#1 #6

PFIHS3
#1 #6

CEMENT SPATULA



CS6
#3

CS5
#3

CS4
#3 #4

CS3
#3

GINGIVAL CORD PACKER



Serrated
GCP113
#3 #4


Non-Serrated
GCP113N
#3 #4

Serrated
GCPAS113
#3 #4

Serrated
GCPR7
#3 #4

Non-Serrated
GCPR7N
#3 #4

IMPRESSION TRAYS DENTULOUS PERFORATED



DENTULOUS NON-PERFORATED



IMPRESSION TRAYS EDENTULOUS PERFORATED

UPPER



ITRLEDPU1

ITRLEDPU2

ITRLEDPU3

ITRLEDPU4

LOWER



ITRLEDPL1

ITRLEDPL2

ITRLEDPL3

ITRLEDPL4

EDENTULOUS NON-PERFORATED

UPPER



ITRLEDNPU1

ITRLEDNPU2

ITRLEDNPU3

ITRLEDNPU4

LOWER



ITRLEDNPL1

ITRLEDNPL2

ITRLEDNPL3

ITRLEDNPL4

IMPRESSION TRAY DENTULOUS PERFORATED PEEDO

UPPER



ITRLDPPU1

ITRLDPPU2

ITRLDPPU3

ITRLDPPU4

LOWER



ITRLDPPL1

ITRLDPPL2

ITRLDPPL3

ITRLDPPL4

IMPRESSION TRAYS HALF

RIGHT



ITRLDPHR1

ITRLDPHR2

ITRLDPHR3

ITRLDPHR4

LEFT



ITRLDPHL1

ITRLDPHL2

ITRLDPHL3

ITRLDPHL4

CROWN REMOVER

TURKEY PATTERN

ANTERIOR AND POSTERIOR

ADJUSTABLE



ITRLDPUA



ITRLDPLA



ITRLDPA



ITRLDNPA



CRT



CRT1



CRT2



CRT3



CRT4



CRT5



CRT6



CRT8



CRT9

CALIPER IWANSON

ARTICULATING PAPER HOLDER



Metal
CLP1



Wax
CLP2



Miller Straight
APF2E



Miller Curved
APF2C

MANUAL PATTERN



CRMP



CRMP1



CRMP2



CRMP3

CROWN REMOVER

AUTOMATIC



CRS

CRP

Crown Remover Tips with CRS & CRP



CR1

CR2

CR3

CR4

CR5

CR6

CR7
(For CRS)

CR8
(For CRP)

CR9

CR10

Basic Kit Includes :

- Crown and Bridge Remover # CRS/CRP
- Single crown remover tip # CR1
- Crown remover tip Blifid # CR2
- Crown & bridge removing tip Sickle # CR3
- Single crown remover Contra angle tip # CR4
- Single crown remover tip Inverted curved # CR5
- Single Crown remover tip straight # CR6
- Bridge remover Flexi-wire set # CR7 (For CRS)
- Bridge remover tip with lock # CR8 (For CRP)
- Tips adjustment Key # CR9
- Open-end adjustment wrench # CR10

Terms of Use

Removal of :

- Permanently cemented crowns and /or bridges
- Temporarily cemented provisional or final restorations
- Maryland bridges
- Orthodontic brackets/bands
- Removable appliances "stuck" at attachments
- Multi-unit restoration "stuck" on working lab model

Precautions :

1. To remove crowns placed over pins or screw type posts, use lowest "strike intensity" to prevent fracture.
2. This instrument could cause damage to the margins of soft gold crowns. If re-use is anticipated, cut a small groove parallel to the lingual margin.
3. Highest strike intensity could cause root fracture; be careful in long span restorations.
4. Only low intensity strike force should be used with teeth having high mobility and little bone support.
5. For highest strike intensity settings, let the individual situation and common sense prevail.

ARTICULATOR



A5222

CROWN HOLDING
FORCEP



CHFRS
12.5cm

CHFRC
12.5cm

PUTTY CUTTER



Right
Handress
PCR

Left
Handress
PCL

PLASTER SHEARS



PCS
21cm

CROWN OPENER



CR4566

TEMPORARY
CROWN REMOVING PLIER



CR4562

Orthodontic

Pliers

Cutters

Black Line Cutters


Accessories



Orthodontic

ORTHODONTICS PLIERS

Every instrument we create becomes, in essence, a critical extension of each practitioner's skill. To achieve this, we handcraft each instrument individually, assuring instruments of the highest quality, durability and performance available to the profession.



Orthodontic

UTILITY PLIERS

Weingart



Short
3000/31



Slim
3000/32



Straight
3000/84



Curved
3000/72
Also available in Wide Tip
4mm Straight # 3000/85
4mm Curved # 3000/73

How

Orthodontic

UTILITY PLIERS



Anterior Band Remover
3000/48



Jhonson Contouring
3000/59



Ball & Socket
3000/4



Band / Crown Crimping
3000/225

Bracket Remover



Curved
3000/82



Straight
3000/83



Short
3000/49



Long
3000/50

Posterior Band Remover



Canine Contouring
3000/310



Separating Module
3000/391



Separating Module Angular
3000/392

Orthodontic

WIRE BENDING PLIERS

Bird Beak



Short
3000/325



Long
3000/326



Tweed O'Brien
3000/97



Three Beak
3000/57

De La Rosa



Wire Contouring
3000/35



Band Contouring
3000/34



Short
3000/51



Long
3000/52

Ribbon/Tweed Arch

Torquing Plier Male/Female Suggested Pair



Female
3000/18



Male
3000/19

Torquing Plier with Key Suggested Pair



Plier
3000/20



Key
3000/21

Orthodontic

WIRE BENDING PLIERS



Nance
3000/80

V Bending
3000/321



Cinch Back
3000/63



Lingual Arch
3000/58

PLIERS WITH TUNGSTEN CARBIDE INSERTS

Band Forming


Band Pinching



Right 3000/42
Left 3000/41



Left 3000/5
Right 3000/6



Weingart TC



Short 3000/31TC
Slim 3000/32TC

LABORATORY PLIERS



Universal 3000/55
Adam's 3000/33
Optical 3000/53



Ligature Tying 3000/210

Ribbon/Tweed Arch

Bird Beak



Short 3000/51TC
Long 3000/52TC
Short TC 3000/325TC
Long TC 3000/326TC

PLIERS WITH TUNGSTEN CARBIDE INSERTS

How TC

4mm



Curved
3000/73TC
Also available in Narrow Tip
2mm Straight # 3000/84
2mm Curved # 3000/72

4mm



Straight
3000/85TC

Posterior Band Remover TC



Long
3000/50TC



Short
3000/49TC



Three Beak Plier with Cutter TC
3000/57TC



Adam's Plier TC
3000/33TC



Crimpable Hook Placement
3000/2



Bracket Removing Curved TC
3000/220

LIGATURE CUTTER TUNGSTEN CARBIDE INSERTS



TUNGSTEN CARBIDE INSERTS



Micro Ligature Cutter
3000/110



Ligature Cutter
3000/708



Ligature Cutter
3000/106



Ligature Cutter
3000/107



Lingual Ligature Cutter 45°
3000/111



Ligature Cutter 15°
3000/109

Orthodontic

DISTAL END CUTTER TUNGSTEN CARBIDE INSERTS



Long Handle
3000/67

Short Handle
3000/67S

With Hole Short Handle
3000/67H


Flush Cut Long Handle
3000/68



Micro Angular
3000/100

Micro
3000/101

HEAVY WIRE CUTTER



Hard Wire Cutter
3000/7

Hard Wire Cutter
3000/8

Hard Wire Cutter
3000/9

Orthodontic

BLACK LINE CUTTER

BLACK LINE

- A performance engineered coating.

REDUCED LIGHT REFLECTION

- The black surface not only prevents reflection of the OT light, but also guarantees optimal contrast in the environment during treatment

TC Material

- Tungsten Carbide material insert in the tip to cut the wire easily.

Distal End Cutter
3200/67B

RIVET JOINT

- Consistent smooth functionality — glides smoothly with no wobble.
- True alignment of cutter tips — rattling at the joint means tips do not come together properly.
- Single part construction — no shims to wear out, seize up or corrode.

Laser Engraving of Applicable Wire Size

- All GDC Ortho cutters & pliers are laser engraved with the applicable wire minimums and maximum



Pin Ligature Cutter
3200/106

Pin Ligature Cutter 15°
3200/109

Pin Ligature Cutter
3200/708

Hard Wire Cutter
3200/7


MATHIEU PLIERS



Smaha NHMS Smaha Curved NHMC Fine NHMF With Ring NHMR

Also available in mathieu smaha with ring # NHMSR

BRACKET HOLDING & PLACING TWEEZER



1006 1008 1132

BUCCAL TUBE HOLDER



1134

LIGATURE INSTRUMENTS



Boley Gauge 4982 Band Pusher 4870 Wire Twister TC 15.0cm WTTC Wire Twister WT WT



Ligature Tucker 672/2 #3 #4 Ligature Directional Instrument 672/1 #3 Utility Instrument 673 #3 Elastic Applicator 674 #3

HEIGHT GAUGE WICK TYPE



018
SS 4876/1



022
SS 4876/2




Adjustable
4875



Band Seater
4871

PLIER STAND
PS



Light Weight Positioner (Aluminium)
2817



Dontrix Gauge
2816



Boon's Gauge SS
2812

For different heights for bracket positioning.
3.5mm, 4mm, 4.5mm & 5mm.
Star points sit inside the brackets face
and measure to the incisal edge.



Boon's Gauge with Lid
2811

Distal Bending Instrument

These instruments for the distal bending of the arches



2815





Instruments MANAGEMENT

Stainless steel trays with
locking instruments in silicon,
resistant to corrosion.
GDC trays are simple to use
and improve the protection
and the storage of your
valuables instruments.
Maximize the efficiency and
performance of washing
and sterilization.

Instruments Management

CASSETTES FOR INSTRUMENTS

cassettes



For 4 Instruments
IMSCA4E



For 5 Instruments
IMSCA5E
IMSCA5PPL



For 6 Instruments
IMSCA6E



For 7 Instruments
IMSCA7E



For 10 Instruments
IMSCA10E
IMSCA10PPL



For 6 Instruments Utility
IMSCA6U

CASSETTES FOR INSTRUMENTS



For Extraction Forceps
IMSCDEF



For Luxatip
IMSCDL



For 20 Instruments
IMSCA20E
IMSCA20PPL




For Surgical Instruments
IMSCDS



For Root Elevator
IMSCDRE


INSTRUMENTS KIT




Extraction Forceps #Premium
EFPWC12




Root Elevators #Premium
REPWC9



Surgical Instruments
SWC10



Perio Surgical
PSIWC12



Endo Surgical
ESIWC13

instruments kit

Instruments Management

INSTRUMENTS KIT



Conservative Instruments #Premium
CIWC19
#4


Basic Implant Surgery Kit
IMPBSWC19



Micro Surgery Kit
MSIWC6



Luxatip New
LWC7



Luxatip
LGWC7

Instruments Management

INSTRUMENTS KIT

Composite Instruments Titanium Coated



PMTWC4
#4



Blue
CIBWC6
#4 #6



Gold
CIGWC6
#4 #6

Composite Instruments Titanium Coated



Filling Instruments
FIWC10
#1 #3 #4



Posterior
CIPWC5
#6



Anterior
CIAWC5
#6

Gracey Curettes



GCWC9
#6



Colour Coded
GCWC7
#7



Colour Coded Mini Five
GCMFWC7
#7

Instruments Management

INSTRUMENTS KIT



Osteotomes Flat Angular
OSTMPUAWC6



Osteotomes Flat Straight
OSTMPUWC6

Also available in concave and convex with cassettes

Osteotomes Concave Angulated
OSTMSHAWC6

Osteotomes Convex Angulated
OSTMSPAWC5



Impaction Kit
IMPACTION



Cleft & Palate Instrument Kit
CLEFT



Trauma Kit
TRAUMA

Instruments Management

INSTRUMENTS KIT



Extraction Forceps #Standard
EFSP12



Extraction Forceps #Premium
EFPP12
EFEP12 (Ergonomic forceps)



Extraction Forceps Peedo #Standard
EFPPSP7



Extraction Forceps Peedo #Premium
EFPPPP7



Root Elevator #Standard
RESP9



Root Elevator #Premium
REPP9



Root Elevator with Warwick James
RESP10



Surgical Instruments
SIP10

OFFER PACKAGE STANDARD
OPS



OFFER PACKAGE PREMIUM
OPP



124.

Instruments Management

INSTRUMENTS KIT



Composite Inst. Gold Flexi Thin
CIGFTP6
#6



Composite Inst. Gold
CIGP6 #6
CIBP6 #6
(Also available in composite
instruments titanium blue)



Composite Inst. Anterior Blue
CIAP5
#6



Composite Inst. Posterior Blue
CIPP5
#6



Conservative Instruments
CIP19
#1 #3 #4



Gracey Curette Premium
GCPP9
#6



Perio Scaler
PSP8
#3 #4 #6



Supra Gingival Scaler
SGSP8

125.

Instruments Management

INSTRUMENTS KIT



Sub Gingival Scaler
SGSP7



Sinus Lift Instruments
IMPSLP8
#6



PKT Carver Colour Coded
PKTCCP5



Impression Trays Dentulous Perforated
ITRLDPP8



Impression Trays Edentulous Perforated
ITRLEDPP8



Impression Trays Dentulous Non-Perforated
ITRLDNP8



Impression Trays Edentulous Non-Perforated
ITRLEDNP8

126.

Instruments Management

INSTRUMENTS KIT



PMT
PMTS4
#3



PMT Premium
PMP4
#4



Filling Instruments
FIS12
#3 #4



RCT Pluggers
RCPS6
#3



RCT Spreaders
RCSS6
#3



Periotomes #5mm
PTS3



Manipal Scalars
MS3
#3



Dental Dam Kit
DDK (ADULT)
DDKP (PEEDO)



Gracey Curettes
GCS9
#3



PKT Carvers
PKTS5
#3

127.

Accessories

INSTRUMENTS KIT



Plier Stand
PS



Extraction Forcep Stand
EFS



Endo Box
AEB



Endo Box Premium
AEBP



Cotton Holder
ACH



Spirit Lamp
ASL



Kidney Tray Autoclavable
KT



Endo Box Round
AEBR



Mouth Props with Chain
MP3

Accessories

X-RAY FILM HANGERS



AFH1



AFH6



AFH10



AFH14

NAPKIN HOLDERS



ANH2



ANH1

NOVELTIES



KCEF



KCRE



KCMMH



KCC



KCIMPT



KCSYR